

# PROTEKT® 600 LIFT

## 33600 | BARIATRIC POWER PATIENT LIFT

The high-performance Protekt® 600 Lift is the perfect blend of form and function for safe and easy patient transfers. This electric-powered, high quality patient transfer system provides an efficient method to safely move patients, preventing patient and caregiver injuries. Designed with the facility in mind, the lift is capable of floor lifts and fits underneath most beds in the lowest position, making it the perfect solution for transferring patients in small rooms with low beds.

HCPCS  
E0635



(Optional Scale  
Item # 30300)

**WARRANTY:**  
Electronics 2 years  
Frame 3 years  
Battery 1 year

### SPECIFICATION FOR ITEM # 33600

Weight Capacity	600 lbs.
Unit Weight	126.5 lbs.
Base Length	47.3"
Base Ext. Width - Closed	26.4"
Base Int. Width - Opened	39.3"
Mast Height	51.3"
Front Caster Height	4.7"
Lifting Height - Lowest	21.5"
Lifting Height - Highest	68.6"
Base Opened Turning Diameter	63.4"
Base Closed Turning Diameter	58.4"

Battery	LP-TC12-BPAK
Charger (Plug)	LP-TP7C-ADP
Charging Cradle (No plug)	LP-TBC2-CRDL
Hand Pendant	LP-TH10-PEN

- **6 point spreader bar allows for multiple sling options, safer patient handling and a wider range of transfers.**
- **Capable of lifting from the floor.**
- **Ergonomically designed foot pedal offers hands free opening of base.**
- **Manual emergency lowering provides safe response to power loss.**
- **4" rear casters with brake & 3" dual front casters.**
- **Clip free sling hookups improve safety and reduce downtime.**
- **600 lb. weight capacity.**
- **Modular design allows easy parts replacement.**
- **Advanced electronics with auxiliary up/down lift controls.**

