

Long Term Care Beds

Each Protopia LTC Bed is unique in its class including the highest quality construction, whisper quiet operation, tool free frame expansion and ease of use.



Protekt Protopia™ Ultra Low LTC Bed

- Lowest position of only 7" helps to mitigate injury from accidental roll outs!
- Integrated tool free length extension (80", 84", 88")
- 8 function hand pendant including auto contour function.
- Heavy duty welded steel construction.
- Rotating pivot assist bars or "Q" bars and head & foot boards included.
- Bumper bar protects walls from damage, includes corner mattress retainers.
- Quick release head and foot boards.
- (8) 3" casters: 4 locking and 4 swivel.
- Weight Capacity: 500 lbs.
- Warranty 5 Years Electrical Components, 15 Years on Frame, 1 Year on all other parts



Protekt Protopia™ EXP Expandable, Ultra Low LTC Bed

- Lowest position of only 7" helps to mitigate injury from accidental roll outs!
- Integrated tool free length extension (80", 84", 88") and width expansion (36", 42", 48")
- 8 function hand pendant including auto contour function.
- Heavy duty welded steel construction.
- Rotating pivot assist bars or "Q" bars and head & foot boards included.
- Bumper bar protects walls from damage, includes corner mattress retainers.
- Quick release head and foot boards.
- (8) 3" casters: 4 locking and 4 swivel.
- Weight Capacity: 550 lbs.
- Warranty 5 Years Electrical Components, 15 Years on Frame, 1 Year on all other parts



Protekt Protopia™ BARI EXP Bariatric "Expandable" LTC Bed

- Lowest position of only 9" helps to mitigate injury from accidental roll outs!
- Integrated tool free length extension (80", 84", 88") and width expansion (36", 42", 48")
- 10 function hand pendant, including auto contour, Trendelenberg and Reverse Trend.
- Heavy duty welded steel construction.
- Rotating pivot assist bars or "Q" bars and head & foot boards included.
- Bumper bar protects walls from damage, includes corner mattress retainers.
- Quick release head and foot boards.
- (8) 3" casters: 4 locking and 4 swivel.
- Weight Capacity: 750 lbs.
- Warranty 5 Years Electrical Components, 15 Years on Frame, 1 Year on all other parts

Pivot Assist Bar



"Q" Bar



Soft Touch Rail



Head & Footboards



Walnut



Cherry



Oak

Mobility, Patient Room Equipment, DME

MOBILITY

Wheelchairs
Transport Chairs
Chairs
Walkers
Sit-to-Stand Lifts
Stand Assist



PATIENT ROOM

Fall Pads (foam and beveled)
Electric Patient Lifts
Patient Lift Slings
I.V. Poles
Overbed Tables



BATH SAFETY

Shower Chairs
Transfer Benches
3 in 1 Commodes
Drop Arm Commodes
Bariatric commodes



ABOUT US

Streamline (verb) - "make an organization or process more efficient and effective by employing faster or simpler working methods"

Our Mission

- Provide innovative, quality, cost-effective healthcare equipment.
- Simplify the process of accessing equipment.
- Provide truly accessible clinical and technical service and support.
- Provide transparent and accurate billing.

Innovative, Quality, Cost-Effective Equipment

Our product offering has been developed using our extensive experience in the design, manufacture and distribution of specialty healthcare equipment. We ensure each product we rent or sell provides true innovation and quality to deliver our customers unequalled value.

Equipment Access

At Streamline Medequip we provide the equipment you need, when you need it, and where you need it. Not all equipment and situations are the same so we give you the option of how to access our products. Sales, rent-to-own, and daily rental: what ever makes sense for your business, Streamline will make it happen.

Service and Support

We provide detailed technical and clinical support including hands-on, in person in-servicing. In addition, we work with our manufacturers to ensure we extend the life of products by backing warranties, providing fast repairs and sending maintenance reminders to extended product life. 24/7, however you chose to contact us (phone, text or e-mail) we will immediately get you the product, information or service you require.

Transparent, Accurate Billing

Transparent online account access allows our customers to see exactly what is being billed and how it is being billed. Review warranty and service information on owned equipment as well as information to make informed equipment decisions.

Download Data

(Online rental review)

CUSTOMER ID	PRODUCT	PO# / PATIENT	START	END	SALE TYPE	STATUS	QUANTITY	RATE	SERIAL#	
AAASNF	Protekt Aire 8000BA Bariatric Air System	FACILITY RENT TO OWN	11/30/2022	2/28/2023	RTO	RTO	1	\$15.00	987656	Details
AAASNF	EC3600 - extriCARE 3600 Pump	SMITH, JANE, 203A	12/22/2022	1/3/2023	DAILY RENTAL	CLOSED	1	\$25.00	123456	Details
AAASNF	EC3600 Drainage canister	SMITH, JANE, 203A	12/22/2022	12/22/2022	SALE	SALE	2	\$25.00		Details
AAASNF	EC3600 Foam Kit	SMITH, JANE, 203A	12/22/2022	12/22/2022	SALE	SALE	6	\$35.00	N/A	Details
AAASNF	EC3600 HydroferaBLUE	SMITH, JANE, 203A	12/29/2022	12/29/2022	SALE	SALE	2	\$25.00		Details
AAASNF	Protekt Aire 6000 LAL/APH System	FACILITY RENT TO OWN	1/2/2023	4/2/2023	RTO	RTO	1	\$8.00	123321	Details
AAASNF	EC3600 - extriCARE 3600 Pump	DOE, JOHN	1/2/2023		DAILY RENTAL	OPEN	1	\$25.00	87678	Details
AAASNF	EC3600 Drainage canister	DOE, JOHN	1/2/2023	1/2/2023	SALE	SALE	2	\$25.00		Details
AAASNF	EC3600 Foam Kit	DOE, JOHN	1/2/2023	1/2/2023	SALE	SALE	6	\$38.00		Details



Streamline Your Access to Specialty Healthcare Equipment



Specialty Mattresses & Cushions



N. P. W. T. (Wound VAC's)



Bariatric Patient Management



Respiratory & Sleep Therapy



Long Term Care Beds



Streamline Advantage Service and Support

Quality Products
Clinical and Technical Support
Personal Service
Rental, Rent-to-Own, Sales
Transparent Rental Review and Billing

www.streamlinemedequip.com

346-445-4017

Specialty Mattress Systems

Streamline MedEquip provides a full continuum of proven Alternating Low-Air-Loss mattress systems. All of our support surfaces provide 8" of air support with low shear/friction patient surfaces and simple operation.



All systems available with foam or air bolsters.


Protekt 4000 - Digital AP/LAL System. Microprocessor control with digital weight in-put for continuously monitored pressure redistribution.

Protekt 6000 - Deluxe Digital AP/LAL System. Advanced digital microprocessor control with adjustable cycle times and fowler boost.

Protekt 8000 - Bariatric Deluxe Digital AP/LAL System. 42" wide surface and 750 lb. Capacity



Pressure Reduction Foam Mattresses



All available with defined perimeter

Protekt pressure redistribution foam mattresses provide three levels of advanced pressure redistribution. Constructed with high spec foam, durable non-skid base and low shear/friction patient surface. Fluid resistant; anti-microbial/anti-bacterial.

Protekt 100 - Low to Moderate risk, single density high specification foam.

Protekt 300 - Multi-tiered die-cut foam in the torso and head sections offers unparalleled pressure redistribution. Sloped heel section.

Protekt 500 with (VGPT)™ - Visco Gel Polymer Technology for rapid heat dissipation and increased breathability. Sloped heel section.

Wheelchair Cushions

Protekt Cushions deliver advanced pressure redistribution and available in a variety of materials and sizes. All covers have durable non-skid bottoms, low shear/friction stretchable patient surface and are fluid resistant; anti-microbial/anti-bacterial.



Protekt® Gel - High density/high resilient foam combined with a gel bladder allows for deep immersion providing excellent pressure redistribution. Dual compartment gel bladder design eliminates gel migration. Safety buckles keeps the cushion in place. Available in all sizes.

Protekt® Air O2 (ROHO Style) - Interconnected air cells provide advanced pressure redistribution. Adjustable 2" or 4" high interconnected air cells offer deep immersion to provide excellent pressure redistribution for tissue protection and treatment. Dual valve compartments address seating challenges such as pelvic obliquities, sacral seating and enhances posture control. Four-way stretch cover enhances immersion and comfort.

Negative Pressure Wound Therapy (Wound VAC's)

Streamline's continuum of NPWT products simplify wound care, improve patient outcomes, assist clinicians with efficiency and reduce overall costs. With a complete line of dressings, accessories and clinical support Streamline delivers exactly what you need to provide advance affordable Negative Pressure Wound Therapy.



extriCARE 3000

Advanced NPWT with LCD touch screen operation, simple settings, on-screen exudate tracking and usage log. Unique intelligent pressure adjustments provide gradual, gentle suction for patient comfort and precise therapy.



extriCARE 3600

Simple, reliable, full function the 3600 is our most affordable VAC. Continuous and intermittent modes of operation, full on screen display provides target and actual pressure for full clinician information.



extriCARE 1000

Small powerful and semi-disposable the 1000 provides many features of larger units for wounds with mild to moderate exudate.

FOAM DRESSING KITS

- 2 step drape and easy to follow instructions (2each)
- Can customize tubing length to tailor to patient needs
- Low profile soft dome port pad
- High tensile strength foam to reduce flaking
- Large, one size fits all kit, (small and medium sizes are available if needed)
- Ruler sticker with each kit to assist with wound Measurement



BANDAGE DRESSINGS

- Anatomically fitted bandage dressings available in various sizes for different wound types
- Non-adhesive contact layer ensures a tight seal and allows negative pressure to reach the wound bed for optimal healing
- Patented flat, perforated drainage tube ensures even distribution of negative pressure across the wound bed
- Bandage Dressing Foam Dressing Connects via tubing to the ExtriCARE canister for effective removal of exudate, maintaining a cleaner, healthier wound site



Mini Oval EC2400-MN-E Small Oval EC2400-SO-E Large Oval EC2400-LO-E Extra Large Oval EC2400-IM-E

ADDITIONAL ACCESSORIES

- Hydrophilic Foam (White Foam)
- Individually Packaged 2 Step Drapes
- No-Sting Skin Prep
- Silver Contact Layer
- Y- Connectors
- Stoma Paste
- Individually Packaged Suction Domes
- Pre-formed Bridges
- Hy-Tape

Bariatric Patient Management

Streamline MedEquip provides a full bariatric continuum to assist in the safe management of bariatric patients.



Protopia EXP

With a 550lb. Capacity and built-in length and width expansion the Protopia EXP accommodates most bariatric patient needs. With a low height of just 7" it can also be used in fall risk management.

Protopia BARI-EXP

When additional capacity is needed the BARI-EXP provides a 750 lb. capacity with a working range of 9" to 31".

Both beds can be provided with an air mattress system or high specification foam mattress as well as half rails or pivot assist bars.



Protekt 8000 Alternating Low-Air-Loss

Deluxe digital control provides 4 cycle times, max inflate and fowler boost. 42" width and 700 lb. capacity accommodates most patient needs.



Protekt 600 Lift and Slings

Patients to 600 lb.



Bariatric Wheelchairs

Streamline provides 22" to 30" bariatric wheelchairs with capacity up to 700 lb.

Wide Reclining (Hi-Back) Wheelchairs

22" & 24" reclining chairs for patients to 350 lb.



Walkers & BSC

Streamline also provides basic bariatric equipment. Regular and drop arm commodes, wheeled walkers and other bariatric DME.

Respiratory & Sleep Therapy



O2 Concentrators

The Alps O2 oxygen 5 liter concentrator is designed for quiet continuous operation. With just one filter maintenance is simple and easy to complete.

Nebulizers

Delivers particles ranging from 0.5 to 5 microns for simple effective therapy.

Floor Suction and Compressors

BiPAP Therapy

The Hypnus BiPAP provides reliable comfortable therapy with four therapy options. BiPAP-S and Auto-S (Spontaneous) BiPAP-T (Timed) BiPAP-ST (Spontaneous Timed)

CPAP Therapy

The APEX CPAP provides basic and auto CPAP with simple reliable operation.



Streamline provides a full line of respiratory and sleep therapy accessories

PAP Masks (Full face, nasal, nasal pillow)
O2 adaptors for CPAP O2 bleed-in
Cabinet and internal filters for most O2 concentrators

Concentrator Maintenance and Quality Assurance

Streamline can provide annual or semi annual maintenance on facility owned concentrators. Service includes external and internal filter replacement, O2 concentration and flow checks.